

## D(-)-炔诺孕酮( $\geq 99\%$ , BioReagent)


产品编号	产品名称	包装
ST1607-50mg	D(-)-炔诺孕酮( $\geq 99\%$ , BioReagent)	50mg
ST1607-250mg	D(-)-炔诺孕酮( $\geq 99\%$ , BioReagent)	250mg
ST1607-1g	D(-)-炔诺孕酮( $\geq 99\%$ , BioReagent)	1g

### 产品简介:

CAS Number	Chemical Formula	Molecular Weight	Purity	Grade
797-63-7	C <sub>21</sub> H <sub>28</sub> O <sub>2</sub>	312.45	$\geq 99\%$	BioReagent

### ➤ 基本信息(General Information):

Name (Chinese)	D(-)-炔诺孕酮
Name (English)	D(-)-Norgestrel
Specifications	$\geq 99\%$
Chemical Formula	C <sub>21</sub> H <sub>28</sub> O <sub>2</sub>
Synonym (Chinese)	-
Synonym (English)	13 $\beta$ -Ethyl-17 $\alpha$ -ethynyl-17 $\beta$ -hydroxygon-4-en-3-one ; 17 $\alpha$ -Ethynyl-18-homo-19-nor-testosterone; 18,19-Dinor-13 $\beta$ -ethyl-17 $\beta$ -hydroxy-4-pregnen-20-yn-3-one; 18,19-Dinor-4-pregnen-20-yn-3-one; Levonorgestrel
Beilstein Registry No.	-
EINECS Number	212-349-8
MDL Number	MFCD00199013
UNSPSC Code	-




### ➤ 产品描述(Description):

Application	D(-)-Norgestrel may be used as a reference standard for the analysis of the analyte in human plasma using liquid chromatography coupled with tandem mass spectrometry (LC-MS/MS).
General description	D(-)-Norgestrel is a hormonal medication and a contraceptive used by women.

### ➤ 性质(Properties):

Assay	$\geq 99\%$
Storage temp.	2-8°C

### ➤ 安全信息(Safety Information):

Hazard Pictogram Codes	
Signal Word	Danger
Hazard Statements	H351-H360D-H362
Precautionary Statements	P203-P260-P263-P264-P270-P308+P316
Personal Protective Equipment	Eyeshields, Gloves, type P3 (EN 143) respirator cartridges
Hazard Codes (Europe)	-
Risk Codes (Europe)	-
Safety Codes (Europe)	-
RIDADR	NONH for all modes of transport
WGK Germany	3

RTECS	-
Flash Point (F)	-
Flash Point (C)	-

**包装清单:**

产品编号	产品名称	包装
ST1607-50mg	D(-)-炔诺孕酮(≥99%, BioReagent)	50mg
ST1607-250mg	D(-)-炔诺孕酮(≥99%, BioReagent)	250mg
ST1607-1g	D(-)-炔诺孕酮(≥99%, BioReagent)	1g
—	说明书	1份

**保存条件:**

4°C保存, 两年有效。

**注意事项:**

- 本产品仅限于专业人员的科学研究用, 不得用于临床诊断或治疗, 不得用于食品或药品, 不得存放于普通住宅内。
- 为了您的安全和健康, 请穿实验服并戴一次性手套操作。

Version 2024.06.27