



碧云天生物技术/Beyotime Biotechnology  
订货热线: 400-1683301 或 800-8283301  
订货 e-mail: order@beyotime.com  
技术咨询: info@beyotime.com  
网址: http://www.beyotime.com

## Trypan blue

产品编号	产品名称	包装
ST798	Trypan blue	1g

### 产品简介:

- Trypan blue也称Direct blue 14, 中文名为台盼蓝。分子式为 $C_{34}H_{24}N_6Na_4O_{14}S_4$ , 分子量为960.81, CAS Number 72-57-1, Cell culture tested, Dye content ~40%。进口分装。
- 本产品为深棕色粉末。
- 常用于细胞死活的检测(cell viability assay)。
- 溶于水, 可以配制成10mg/ml的水溶液, 水溶液呈深蓝色。

### 包装清单:

产品编号	产品名称	包装
ST798	Trypan blue	1g
—	说明书	1份

### 保存条件:

室温保存。

### 注意事项:

- 本产品对人体有毒, 操作时请特别小心, 并注意有效防护以避免直接接触人体或吸入体内。
- 本产品仅限于专业人员的科学研究用, 不得用于临床诊断或治疗, 不得用于食品或药品, 不得存放于普通住宅内。
- 为了您的安全和健康, 请穿实验服并戴一次性手套操作。

### 使用本产品的文献:

1. Liu Z, Jiang M, Kang T, Miao D, Gu G, Song Q, Yao L, Hu Q, Tu Y, Pang Z, Chen H, Jiang X, Gao X, Chen J. . Lactoferrin-modified PEG-co-PCL nanoparticles for enhanced brain delivery of NAP peptide following intranasal administration. *Biomaterials*. 2013 May;34(15):3870-81.
2. Huang P, Pandoli O, Wang XS, Wang Z, Li ZM, Zhang CL, Chen F, Lin J, DX Cui, Chen XY. . Chiral guanosine 5'-monophosphate-capped gold nanoflowers: Controllable synthesis, characterization, surface-enhanced Raman scattering activity, cellular imaging and photothermal therapy. *NANO RESEARCH*. 2012 Sep;5(9):630-9.
3. Zhou ZH, Song JW, Li W, Liu X, Cao L, Wan LM, Tan YX, Ji SP, Liang YM, Gong F. . The acid-sensing ion channel, ASIC2, promotes invasion and metastasis of colorectal cancer under acidosis by activating the calcineurin/NFAT1 axis. *J EXP CLIN CANC RES*. 2017 Sep 19;36(1):130.
4. Liu F, Cheng Z, Li X, Li Y, Zhang H, Li J, Liu F, Xu H, Li F. . A Novel Pak1/ATF2/miR-132 Signaling Axis Is Involved in the Hematogenous Metastasis of Gastric Cancer Cells. *MOL THER-NUCL ACIDS*. 2017 Sep 15;8:370-382.
5. Zhang C, Lei JL, Zhang H, Xia YZ, Yu P, Yang L, Kong LY. . Calyxin Y sensitizes cisplatin-sensitive and resistant hepatocellular carcinoma cells to cisplatin through apoptotic and autophagic cell death via SCF  $\beta$ TrCP-mediated eEF2K degradation. *ONCOTARGET*. 2017 Aug 3;8(41):70595-70616.
6. Chang Y, Yang ST, Liu JH, Dong E, Wang Y, Cao A, Liu Y, Wang H. . In vitro toxicity evaluation of graphene oxide on A549 cells. *Toxicol Lett*. 2011 Feb 5;200(3):201-10.
7. Wu L, Chen SX, Jiang XG, Fu XL, Shen YJ, Cao JP. . Separation and purification of *Toxoplasma gondii* tachyzoites from in vitro and in vivo culture systems. *Exp Parasitol*. 2012 Jan;130(1):91-4.
8. Shen XC, Tao L, Li WK, Zhang YY, Luo H, Xia YY. . Evidence-based antioxidant activity of the essential oil from *Fructus A. zerumbet* on cultured human umbilical vein endothelial cells' injury induced by ox-LDL. *BMC COMPLEMENT ALTERN M*. 2012 Oct 7;12:174.
9. Xia H, Gu G, Hu Q, Liu Z, Jiang M, Kang T, Miao D, Song Q, Yao L, Tu Y, Chen H, Gao X, Chen J. . Activatable cell penetrating peptide-conjugated nanoparticles with enhanced permeability for site-specific targeting delivery of anticancer drug. *BIOCONJUGATE CHEM*. 2013 Mar 20;24(3):419-30.
10. Jie Q, Tang Y, Deng Y, Li Y, Shi Y, Gao C, Xing M, Wang D, Liu L, Yao P. . Bilirubin participates in protecting of heme oxygenase-1 induction by quercetin against ethanol hepatotoxicity in cultured rat hepatocytes. *Alcohol*. 2013 Mar;47(2):141-8.
11. Shen J, Wang J, Shao YR, He DK, Zhang L, Nadeem L, Xu G. . Adenovirus-delivered angiotensin-1 treatment for phosgene-induced acute lung injury. *Inhal Toxicol*. 2013 Apr;25(5):272-9.
12. Liu Z, Gao X, Kang T, Jiang M, Miao D, Gu G, Hu Q, Song Q, Yao L, Tu Y, Chen H, Jiang X, Chen J. . B6 Peptide-Modified PEG-PLA Nanoparticles for Enhanced Brain Delivery of Neuroprotective Peptide. *BIOCONJUGATE CHEM*. 2013 Jun 19;24(6):997-1007.
13. He Y, Zhou JW, Xu L, Gong MJ, He TC, Bi Y. . Comparison of proliferation and differentiation potential between mouse primary hepatocytes and

embryonic hepatic progenitor cells in vitro. *Int J Mol Med*. 2013 Aug;32(2):476-84.

14. Yang YH, Chen K, Li B, Chen JW, Zheng XF, Wang YR, Jiang SD, Jiang LS. . Estradiol inhibits osteoblast apoptosis via promotion of autophagy through the ER-ERK-mTOR pathway. *Apoptosis*. 2013 Nov;18(11):1363-75.
15. Chen Y, Wang G, Liu Z, Wang S, Wang Y. . Glucocorticoids regulate the proliferation of T cells via miRNA-155 in septic shock. *Exp Ther Med*. 2016 Dec;12(6):3723-3728.
16. Huang Y, Wang XX, Sun DD, Zhang ZX, Yang WW, Shao T, Han H, Zhang EF, Pu ZS, Hou ZX, Dong HL, Xiong LZ, Hou LC. . Sub-anesthesia Dose of Isoflurane in 60% Oxygen Reduces Inflammatory Responses in Experimental Sepsis Models. *CHINESE MED J-PEKING*. 2017 Apr 5;130(7):840-853.
17. Li J, Xu Z, Chupalov A, Marchisio MA . Anti-CRISPR-based biosensors in the yeast *S. cerevisiae*. *J Biol Eng*. 2018 Aug 1 12:11.
18. Zhao X, Jin Y, Yang L, Hou Z, Liu Y, Sun T, Pei J, Li J, Yao C, Wang X, Chen G . Promotion of SIRT1 protein degradation and lower SIRT1 gene expression via reactive oxygen species is involved in Sb-induced apoptosis in BEAS-2b cells. *Toxicol Lett*. 2018 Oct 15 296:73-81.
19. Meng G, Fei Z, Fang M, Li B, Chen A, Xu C, Xia M, Yu D, Wei J . Fludarabine as an Adjuvant Improves Newcastle Disease Virus-Mediated Antitumor Immunity in Hepatocellular Carcinoma. *Mol Ther Oncolytics*. 2019 Mar 27 13:22-34.
20. Zhang L, Wang H, Zhou X, Mao L, Ding K, Hu Z . Role of mitochondrial calcium uniporter-mediated Ca<sup>2+</sup> and iron accumulation in traumatic brain injury. *J Cell Mol Med*. 2019 Apr 23(4):2995-3009.

Version 2024.03.12