



碧云天网站



微信公众号

碧云天生物技术/Beyotime Biotechnology

订货热线: 400-168-3301或800-8283301

订货e-mail: order@beyotime.com

技术咨询: info@beyotime.com

网址: http://www.beyotime.com

Pifithrin- α (p53抑制剂)

产品编号	产品名称	包装
S1816-5mg	Pifithrin- α (p53抑制剂)	5mg
S1816-25mg	Pifithrin- α (p53抑制剂)	25mg

产品简介:

- 本Pifithrin- α 为Pifithrin- α ·hydrobromide, 简称Pft α , 是一种常用的p53抑制剂。Pifithrin- α 能可逆性地抑制p53依赖的基因转录激活, 例如cyclin G, p21/waf1和mdm2; 也可以抑制p53介导的细胞凋亡。
- 本产品可以用于细胞培养, 也可以用于动物实验。
- Pifithrin- α ·hydrobromide的分子量为367.30, 分子式为C₁₆H₁₈N₂O₅·HBr, CAS Number: 63208-82-2。本产品为进口分装, 纯度大于98%。
- Pifithrin- α 可溶于DMSO或甲醇。

包装清单:

产品编号	产品名称	包装
S1816-5mg	Pifithrin- α (p53抑制剂)	5mg
S1816-25mg	Pifithrin- α (p53抑制剂)	25mg
—	说明书	1份

保存条件:

-20°C避光保存。

注意事项:

- 本产品仅限于专业人员的科学研究用, 不得用于临床诊断或治疗, 不得用于食品或药品, 不得存放于普通住宅内。
- 为了您的安全和健康, 请穿实验服并戴一次性手套操作。

使用说明:

1. Pifithrin- α 常见使用浓度范围为5-20 μ M。实际的最佳工作浓度请参考相关文献, 或根据实验目的, 以及所培养的特定细胞和组织, 通过实验进行摸索和优化。

使用本产品的文献:

1. Han YT, Chen XH, Xie J, Zhan SM, Wang CB, Wang LX. Purple sweet potato pigments scavenge ROS, reduce p53 and modulate Bcl-2/Bax to inhibit irradiation-induced apoptosis in murine thymocytes. Cell Physiol Biochem. 2011;28(5):865-72.
2. Wang Y, Liu Q, Liu Z, Li B, Sun Z, Zhou H, Zhang X, Gong Y, Shao C. Berberine, a genotoxic alkaloid, induces ATM-Chk1 mediated G2 arrest in prostate cancer cells. Mutat Res. 2012 Jun 1;734(1-2):20-9.
3. Yang M, Yuan F, Li P, Chen Z, Chen A, Li S, Hu C. Interferon regulatory factor 4 binding protein is a novel p53 target gene and suppresses cisplatin-induced apoptosis of breast cancer cells. Mol Cancer. 2012 Aug 13;11:54.
4. Qiao H, Chen X, Xu L, Wang J, Zhao G, Hou Y, Ge HM, Tan RX, Li E. Antitumor effects of naturally occurring oligomeric resveratrol derivatives. FASEB J. 2013 Nov;27(11):4561-71.
5. Wang J, Niu Z, Shi Y, Gao C, Wang X, Han J, Li J, Gao Z, Zhu X, Song X, Qin Z, Wang H. Bcl-3, induced by Tax and HTLV-1, inhibits NF- κ B activation and promotes autophagy. Cell Signal. 2013 Dec;25(12):2797-804.
6. Guo H, Liu Z, Xu B, Hu H, Wei Z, Liu Q, Zhang X, Ding X, Wang Y, Zhao M, Gong Y, Shao C. Chemokine receptor CXCR2 is transactivated by p53 and induces p38-mediated cellular senescence in response to DNA damage. Aging Cell. 2013 Dec;12(6):1110-21.
7. Xu B, Wang W, Guo H, Sun Z, Wei Z, Zhang X, Liu Z, Tischfield JA, Gong Y, Shao C. Oxidative stress preferentially induces a subtype of micronuclei and mediates the genomic instability caused by p53 dysfunction. Mutat Res Fundam Mol Mech Mutagen. 2014 Dec 1;770:1-8.
8. Song X, Li L, Shi Q, Lehmler HJ, Fu J, Su C, Xia X, Song E, Song Y. Polychlorinated Biphenyl Quinone Metabolite Promotes p53-Dependent DNA Damage Checkpoint Activation, S-Phase Cycle Arrest and Extrinsic Apoptosis in Human Liver Hepatocellular Carcinoma HepG2 Cells. Chem Res Toxicol. 2015 Nov 16;28(11):2160-9.
9. Liu YP, Dong FX, Chai X, Zhu S, Zhang BL, Gao DS. Role of Autophagy in Capsaicin-Induced Apoptosis in U251 Glioma Cells. Cell Mol Neurobiol. 2016 Jul;36(5):737-43.
10. Ren J, Li XR, Liu PC, Cai MJ, Liu W, Wang JX, Zhao XF. G-protein α q participates in the steroid hormone 20-hydroxyecdysone nongenomic signal transduction. J Steroid Biochem Mol Biol. 2014 Oct;144 Pt B:313-23.
11. Chen J, Jiang L, Lan K, Chen X. Celecoxib Inhibits the Lytic Activation of Kaposi's Sarcoma-Associated Herpesvirus through Down-Regulation of RTA Expression by Inhibiting the Activation of p38 MAPK. Viruses. 2015 May 5;7(5):2268-87.

Version 2016.12.12