

# Honokiol up-regulates prostacyclin synthase protein expression and inhibits endothelial cell apoptosis

Xiaoxue Zhang<sup>a</sup>, Shizhong Chen<sup>b</sup>, Yinye Wang<sup>a,\*</sup>

<sup>a</sup> Department of Pharmacology, School of Pharmaceutical Sciences, Peking University, Beijing 100083, China

<sup>b</sup> Department of Natural Medicinal Chemistry, Peking University, Beijing 100083, China

Received 4 April 2006; received in revised form 26 September 2006; accepted 27 September 2006

Available online 10 October 2006

## Abstract

Honokiol is a bioactive compound extracted from the Chinese medicinal herb *Magnolia officinalis*. We recently demonstrated that honokiol inhibited arterial thrombosis through stimulation of prostacyclin (PGI<sub>2</sub>) generation and endothelial cell protection. The current study is designed to investigate its mechanism of stimulation of PGI<sub>2</sub> generation and cell protection. 6-keto-PGF1 $\alpha$ , the stable metabolite of PGI<sub>2</sub>, in the media of rat aortic endothelial cells was measured with radioimmunoassay kits. Indomethacin, an inhibitor of cyclooxygenase (COX) and tranlycypromine, a prostacyclin synthase inhibitor were used to ascertain the target enzyme affected by honokiol. Prostacyclin synthase protein levels in endothelial cells were determined by Western blot analysis using an anti-PGI<sub>2</sub> synthase rabbit polyclonal antibody. Flow cytometry was used to quantify the apoptotic cells and spectrophotometry was used to test the caspase-3 activity. Honokiol (0.376–37.6  $\mu$ M) increased the level of 6-keto-PGF1 $\alpha$  in the media of normal endothelial cells. It counteracted the inhibitory effect of tranlycypromine on the PGI<sub>2</sub> generation, but did not influence the effect of indomethacin; evidently, honokiol up-regulated the expression of prostacyclin synthase in the endothelial cells. These effects showed perfect concentration-dependent behavior. In addition, at lower concentration (0.376–3.76  $\mu$ M), honokiol significantly decreased the percentage of apoptotic endothelial cells induced by oxidized low-density lipoprotein (ox-LDL) and significantly lowered the activity of caspase-3 stimulated by ox-LDL. A high dose of honokiol (37.6  $\mu$ M), however, failed to influence either of them. In conclusion, honokiol augments PGI<sub>2</sub> generation in normal endothelial cells; its effect on PGI<sub>2</sub> generation attributes to up-regulation of prostacyclin synthase expression; its cell protection may be correlated with its inhibition on apoptosis of endothelial cells. These findings have partly revealed the molecular mechanism of honokiol on inhibiting arterial thrombosis.

© 2006 Elsevier B.V. All rights reserved.

**Keywords:** Honokiol; Thrombosis; Rat aortic endothelial cells; Prostacyclin (PGI<sub>2</sub>); Prostacyclin synthase; Apoptosis

## 1. Introduction

Thrombosis is a major cause of cardiovascular and cerebrovascular diseases. Endothelial cells play important roles in controlling thrombosis by synthesis and release of endogenous active substances such as prostacyclin (PGI<sub>2</sub>), nitric oxide (NO), tissue-plasminogen, etc. It is these substances that are directly involved in regulating blood flow and vascular stress. Structural and functional integrity of endothelial cells are important in the maintenance of vessel walls and circulatory function. Endothelial dysfunction facilitates a number of unfavorable effects on vasculature, including platelet activation on injured endotheli-

um (Cines et al., 1998), which then leads to an imbalance in endothelium-derived factors and, thereafter, induction of the formation and development of thrombi (Golino et al., 1989; Kishi and Numano, 1989; Watanabe et al., 1997).

Honokiol is a main biphenyl neolignan isolated from traditional Chinese medicine Hou p'u, cortex of *Magnolia officinalis* (*Magnoliaceae*), which has a variety of pharmacological effects such as anti-oxidation (Lo et al., 1994; Chiu et al., 1997), anti-platelet aggregation (Teng et al., 1988), anxiolytic effect (Watanabe et al., 1983), etc. Recently, we demonstrated that honokiol inhibited arterial thrombosis through endothelial cell protection and stimulation of PGI<sub>2</sub> generation (Hu et al., 2005), but the detailed mechanisms of these effects were not clear yet. In current study we explored its action mechanism on the generation of PGI<sub>2</sub> and on cell protection, further explaining honokiol's anti-thrombotic effects.

\* Corresponding author. Tel.: +86 10 8280 2652; fax: +86 10 6201 5584.

E-mail address: [wangyinye@bjmu.edu.cn](mailto:wangyinye@bjmu.edu.cn) (Y. Wang).

## 2. Materials and methods

### 2.1. Materials

Honokiol injection was prepared by the Department of Natural Medicinal Chemistry, Peking University. The water-soluble preparation of honokiol contains polyvinyl and other auxiliary materials. RPMI Medium 1640 was the product of Invitrogen Corporation, USA; NBCS (Newborn Calf Serum) was the product of PAA Laboratoies GmbH; and ox-LDL was the product of Beijing Union Sanyou Science & Technology Development Corporation Ltd, China. Rnase, indomethacin, arachidonic acid and tranlycypromine were all products of Sigma, USA; 6-keto-PGF $_{1\alpha}$  immunoassay kits were products of Radioimmunity Institute of People's Liberation Army General Hospital, Beijing; Cell lysis buffer was the product of Beyotime Biotechnology, China; and the vehicle control was normal saline plus definite proportion of auxiliary materials.

### 2.2. Cell culture

Male Sprague–Dawley rats weighing 180–200 g were anesthetized with an overdose of sodium pentobarbital and the abdominal aortas of rats were rapidly removed and collected in 1640 medium. Surrounding fat and connective tissue were cleaned off, and then aortas were opened longitudinally. Cells from the endothelium were harvested by gently scraping the luminal surface using vertical ophthalmic forceps. The cells were then seeded in T75 polystyrene flasks pre-coated with sterile murine tail collagen in the medium 1640 containing 20% Newborn Calf Serum, 100 U/ml of penicillin and 100  $\mu$ g/ml of streptomycin, and cultured at 37 °C in a humidified 5% CO $_2$  atmosphere. Cells were allowed to grow undisturbed for 4 days and thereafter the media were changed every 2 days for a total culturing period of 8 to 10 days. The endothelial cells were identified according to the previous method (Jaffe et al., 1973) and were subcultured as described by others (Centra et al., 1992). Experiments were performed using confluent cultures of passage 5 (Sibel et al., 2001).

### 2.3. Measurement of 6-keto-PGF $_{1\alpha}$ in media of normal rat aortic endothelial cells

The rat aortic endothelial cells were seeded into 96-well plates with  $1 \times 10^5$  cells/ml and cultured until confluence. The media were removed and cells were washed twice with PBS, and pre-incubated for 30 min with the medium 1640 containing 2% Newborn Calf Serum. The cultures were then incubated for 24 h with vehicle, 0.376, 3.76 or 37.6  $\mu$ M of honokiol respectively. Culture media were saved for the determination of 6-keto-PGF $_{1\alpha}$ , which was determined with immunoassay kits.

### 2.4. Intervention experiments of inhibitors of cyclooxygenase (COX) and prostacyclin synthase

The rat aortic endothelial cells were seeded into 96-well plates with  $1 \times 10^5$  cells/ml and cultured until confluence. The

media were removed and cells were washed twice with PBS, and then pre-incubated with the medium 1640 containing 2% Newborn Calf Serum for 30 min. After pre-treated with vehicle, 0.376, 3.76 or 37.6  $\mu$ M of honokiol at 37 °C for 30 min, the cultures were incubated with indomethacin (3  $\mu$ M) for 24 h (Ottino and Duncan, 1997) to evaluate the influence of honokiol on COX. In the other set of samples, after treated with vehicle, 0.376, 3.76 or 37.6  $\mu$ M of honokiol at 37 °C for 30 min, the cultures were incubated with the prostacyclin synthase inhibitor tranlycypromine (1 mM) for 24 h to evaluate the influence of honokiol on prostacyclin synthase (Debbie and Mary, 1995). The vehicle controls and inhibitor controls were set up to observe their influences. Arachidonic acid (30  $\mu$ M) was added to all cultures to supply sufficient substrate. Culture fluids were saved for the determination of PGI $_2$ , which is the final product of COX or prostacyclin synthase. 6-keto-PGF $_{1\alpha}$  in cultured fluid was determined with immunoassay kits.

### 2.5. Western blot analysis of prostacyclin synthase

The rat aortic endothelial cells were seeded into 6-well plates with  $2 \times 10^5$  cells/ml and cultured until confluence. Endothelial cells were treated with tumor necrosis factor- $\alpha$  (TNF $\alpha$ ) (10 ng/ml), vehicle, 0.376, 3.76 or 37.6  $\mu$ M of honokiol for 24 h, respectively. After removal of the media, cells were washed twice with ice-cold PBS, then lysed using cell lysis buffer [20 mM Tris PH7.5, 150 mM NaCl, 1% Triton X-100, 2.5 mM sodium pyrophosphate, 1 mM EDTA, 1% Na $_3$ VO $_4$ , 0.5  $\mu$ g/ml leupeptin, 1 mM phenylmethanesulfonyl fluoride (PMSF)]. The lysates were collected by scraping from the plates, and then centrifuged at 10,000  $\times g$  at 4 °C for 5 min. Western blot was performed according to the procedure of Towbin et al. (1979). Briefly, proteins (30  $\mu$ g) were loaded on a 12% of SDS-polyacrylamide gel for electrophoresis, then transferred onto nitrocellulose transfer membranes (Osmonics, USA) at 0.8 mA/cm $^2$  for 2 h. Membranes were blocked at room temperature for 1 h with blocking solution [5% skimmed milk in Tris-buffered solution plus Tween-20 (TBST): 50 mM Tris-HCl, 150 mM NaCl, pH=7.5, 0.1% v/v Tween 20]. Membranes were then incubated overnight at 4 °C with anti-PGI $_2$  synthase rabbit polyclonal antibody at 1:150 (Santa Cruz Biotechnology, USA) in blocking solution. After two 10 min washing in TBST, membranes were incubated for 1 h at room temperature with a horseradish peroxidase-conjugated anti-rabbit secondary antibody at 1:500 in blocking solution. Detection was performed by enhanced chemiluminescence (ECL) using a Western Blotting Luminol Reagent (Santa Cruz Biotechnology, USA) according to the manufacturer's instructions. Bands were then quantified by scanning densitometry (THERMAL IMAGING SYSTEM FTI-500, Pharmacia Biotech). Protein concentrations were determined using the Pierce Micro BCA protein assay system (Pierce, Rockford, Illinois, USA).  $\beta$ -actin of rat aortic endothelial cells was used as a housekeeping protein, and was determined following the same procedure as mentioned above, using a specific anti-actin mouse monoclonal antibody (Sigma-Aldrich, Madrid, Spain) at 1:1000 in blocking solution and a horseradish peroxidase-




Fig. 1. Effect of honokiol on the level of 6-keto-PGF $\alpha$  in media of normal rat aortic endothelial cells. Rat aortic endothelial cells of passage 5 were seeded with  $1 \times 10^5$  cells/ml in 96-well plates and grown to confluence, then treated with vehicle and honokiol for 24 h, respectively. 6-keto-PGF $\alpha$  in culture fluids was determined with immunoassay kits. Each datum was from 6 repeated tests, and expressed as mean  $\pm$  S.D. \* $P < 0.05$ ; \*\* $P < 0.01$ .

conjugated anti-mouse secondary antibody at 1:500 in blocking solution.

## 2.6. Flow cytometry analysis

The rat aortic endothelial cells were seeded into 60 mm culture dishes with  $5 \times 10^5$  cells/ml and cultured until confluence. The cells were washed twice with PBS and pre-incubated for 30 min with the medium 1640 containing 2% Newborn Calf Serum. The cells were then treated with vehicle, 0.376, 3.76 or 37.6  $\mu$ M of honokiol at 37  $^{\circ}$ C for 30 min, respectively, and incubated with oxidized low-density lipoprotein (ox-LDL) (final concentration was 100  $\mu$ g/ml) for 24 h to induce cell apoptosis. The degree of apoptotic cells was analyzed by flow cytometry as previously described (Chaouchi et al., 1994). Shortly, pelleted cells were suspended and fixed in




Fig. 2. Influence of honokiol on COX inhibited by indomethacin. Rat aortic endothelial cells passage 5 were seeded at  $1 \times 10^5$  cells/ml in 96-well plates and grown to confluence, then were exposed to indomethacin (3  $\mu$ M) for 24 h after treated with vehicle or honokiol for 30 min. 6-keto-PGF $\alpha$  in culture fluids was determined with immunoassay kits. Each datum was from 6 repeated tests, and expressed as mean  $\pm$  S.D. \*\* $P < 0.01$ . "IND" represents indomethacin.




Fig. 3. Influence of honokiol on prostacyclin synthase inhibited by tranylcypromine. Rat aortic endothelial cells of passage 5 were seeded at  $1 \times 10^5$  cells/ml in 96-well plates and grown to confluence, then were exposed to tranylcypromine (1 mM) for 24 h after treated with vehicle or honokiol for 30 min. 6-keto-PGF $\alpha$  in culture fluids was determined with immunoassay kits. Each datum was from 6 repeated tests, and expressed as mean  $\pm$  S.D. \* $P < 0.05$ , \*\* $P < 0.01$ . "TRAN" represents tranylcypromine.

70% cold ethanol overnight at  $-20$   $^{\circ}$ C after twice washing with PBS. Then the cells were re-suspended in PBS containing RNase A (200  $\mu$ g/ml), and incubated at 37  $^{\circ}$ C for 30 min. Propidium iodide PI was added to the cell suspensions to give final concentrations of 100  $\mu$ g/ml. The fluorescence intensity of




Fig. 4. Effect of honokiol on prostacyclin synthase protein expression in rat aortic endothelial cells. Rat aortic endothelial cells of passage 5 were seeded at  $2 \times 10^5$  cells/ml in 6-well plates and grown to confluence, then were exposed to vehicle, TNF $\alpha$  (10 ng/ml) or honokiol for 24 h. Then Western blot analysis of cellular proteins from different treated cells was performed. A) The bands of prostacyclin synthase and  $\beta$ -actin were identified by Western blot analysis using an anti-prostacyclin synthase rabbit polyclonal antibody and anti-actin mouse monoclonal antibody, in which  $\beta$ -actin was used as a loading control. Shown was one Western blot result of 3 separate experiments; B) The statistical analysis results of ratio of prostacyclin synthase compared to  $\beta$ -actin. Results were depicted as mean  $\pm$  S.D. \* $P < 0.05$ , \*\* $P < 0.01$ .

PI was analyzed with a FACScalibur flow cytometer and CellQuest software.

### 2.7. Caspase-3 activity assay

The rat aortic endothelial cells were seeded into 24-well plates with  $2 \times 10^6$  cells/ml and cultured until confluence, and the apoptosis was induced as mentioned above. At the same time, honokiol controls were set up in which cells were treated only with honokiol (0.376, 3.76, and 37.6  $\mu\text{M}$ ) without ox-LDL. After induced apoptosis the activity of caspase-3 in cells was measured with spectrophotometry at 405 nm (MICROPLATE READER MODEL 550, BIO-RAD) using a commercially available kit (Calbiochem) according to the manufacturer's instruction, which detects chromophore *p*-nitroanilide after cleavage from the labeled substrate *N*-acetyl-Asp-Glu-Val-Asp-*p*-nitroanilide.

### 2.8. Statistical analysis

The data are expressed as mean  $\pm$  S.D. Statistical evaluation was performed using Dunnett *t*-test to compare the differences between treated groups and control groups.  $P < 0.05$  was considered to be statistically significant.

## 3. Results

### 3.1. Influence of honokiol on $\text{PGI}_2$ generation in normal rat aortic endothelial cells

Normal endothelial cells treated with 0.376–37.6  $\mu\text{M}$  of honokiol produced a significant concentration-dependent increase in the level of 6-keto- $\text{PGF}_1\alpha$  in media compared with vehicle-treated cells (Fig. 1), it was demonstrated that honokiol augmented the generation of  $\text{PGI}_2$  besides its protective effect on endothelial cells thereby increasing  $\text{PGI}_2$  generation as our previous presupposition (Hu et al., 2005).

### 3.2. Influence of honokiol on COX inhibited by indomethacin

In the presence of COX inhibitor the concentration of 6-keto- $\text{PGF}_1\alpha$  in the media was measured to evaluate the effect of honokiol on COX. The treatment with indomethacin plus vehicle resulted in a significant decrease on the concentration of 6-keto- $\text{PGF}_1\alpha$  in the media. The treatment with indomethacin plus 0.376, 3.76 or 37.6  $\mu\text{M}$  of honokiol failed to change the effect of indomethacin (Fig. 2). The result demonstrated that honokiol did not influence COX in the presence of the inhibitor. In other words, the influence of honokiol on  $\text{PGI}_2$  is irrelative with COX.

### 3.3. Influence of honokiol on prostacyclin synthase inhibited by tranlycypromine

Determining the concentration of 6-keto- $\text{PGF}_1\alpha$  in endothelial cells media was used to assess the influence of honokiol on prostacyclin synthase. The treatment with tranlycypromine,




Fig. 5. Effect of honokiol on rat aortic endothelial cells apoptosis induced by ox-LDL. Rat aortic endothelial cells of passage 5 were seeded at  $5 \times 10^5$  cells/ml in 60 mm culture dishes and grown to confluence. After pre-incubated with vehicle or honokiol for 30 min, cultures were exposed to ox-LDL for 24 h to induce apoptosis; the cells were then collected to determine apoptotic percentages. Apoptotic cells were quantified by FACS analysis of cell nuclei after staining with propidium iodide. Cells containing amounts of DNA found in sub-G1 peak are in region M1, and the ratio of M1 and total region represented the percentage of apoptotic cells. A) Shown was one representative of apoptosis determination from 3 repeated separate experiments. B) The statistical analysis results of percentage of apoptotic cells. Results were depicted as mean  $\pm$  S.D. \*\* $P < 0.01$ .

a prostacyclin synthase inhibitor, resulted in a significant decrease on the concentration of 6-keto-PGF<sub>1</sub>α compared with vehicle. 0.376–37.6 μM of honokiol reversed the influence of prostacyclin synthase inhibitor on 6-keto-PGF<sub>1</sub>α with a concentration-dependent behavior (Fig. 3). The results indicate that prostacyclin synthase was the target enzyme affected by honokiol.

#### 3.4. Effect of honokiol on prostacyclin synthase protein expression

We investigated the effect of honokiol on prostacyclin synthase expression with western blot analysis. TNFα was used as a positive control according to other study (Romualdo et al., 1997). The treatment with honokiol or TNFα increased the prostacyclin synthase protein levels compared with vehicle, and the effect of honokiol exhibited a concentration dependent manner (Fig. 4). This result manifested that honokiol stimulates prostacyclin synthase expression in rat aortic endothelial cells; it may be a promoter of prostacyclin synthase expression.

#### 3.5. Effect of honokiol on apoptosis of rat aortic endothelial cells

Flow cytometry was used to quantify the amount of apoptotic cells. The amount of bound dye correlates with the DNA content in a given cell counts, and DNA fragmentation in apoptotic cells translates into fluorescence intensity lower than that of G1 cells (sub-G1 peak) corresponding to apoptotic cells. The percentage of apoptosis for 0.376 and 3.76 μM of honokiol treated cells were significantly lowered compared with vehicle treated cells, but treatment with 37.6 μM of honokiol failed to influence the percentage of apoptosis (Fig. 5), indicating that at

appropriate concentrations (0.376–3.76 μM) honokiol inhibits endothelial cell apoptosis, which may be one of its mechanisms for cell protection.

#### 3.6. Effect of honokiol on caspase-3 activity of apoptotic rat aortic endothelial cells

Caspase-3 activity induced by ox-LDL in the presence or absence of honokiol was analyzed. As shown in Fig. 6, honokiol (0.376–3.76 μM) significantly inhibited the caspase-3 activity induced by ox-LDL, but a high dose of honokiol (37.6 μM) failed to show effect. However, honokiol did not influence the caspase-3 activity of normal rat aortic endothelial cells, which was consistent with our previous finding that honokiol caused no damage to endothelial cells (Hu et al., 2005).

## 4. Discussion

Honokiol demonstrated platelet aggregation inhibition by inhibiting thromboxane A<sub>2</sub> (TXA<sub>2</sub>) formation in platelets (Teng et al., 1988). Recently it was found to prevent cerebral injury caused by middle cerebral artery occlusion and cerebral ischemia reperfusion injury, the mechanism for these effects maybe its antioxidant activity (Liou et al., 2003; Chen et al., 2003). More recently it was discovered that honokiol significantly inhibited arterial thrombosis induced by endothelium injury in rats, which may be attributed to its protective effect on endothelial cells and its stimulation of PGI<sub>2</sub> generation (Hu et al., 2005). Therefore, it is presupposed that its protective effect on injured rat aortic endothelial cells may be one of the reasons for the increase in the PGI<sub>2</sub> generation. The viability of cells treated with 37.6 μM of honokiol is not highest, but the PGI<sub>2</sub> level in media is the highest at this concentration (Hu et al., 2005), which leads to the conclusion that there should be other mechanisms than cells protection for PGI<sub>2</sub> level increase. In this study, honokiol (0.376–37.6 μM) significantly and concentration-dependently increased PGI<sub>2</sub> level in normal rat aortic endothelial cells media, which further suggested that honokiol increased PGI<sub>2</sub> production through pathways other than the protection of endothelial cells.

COX and prostacyclin synthase are the two rate-limiting enzymes of PGI<sub>2</sub> synthesis in endothelial cells (Feng et al., 1993; Spisni et al., 1995). Therefore, COX inhibitor indomethacin or prostacyclin synthase inhibitor tranlylcypromine were used to explore the effects of honokiol on the two enzymes. The results showed that honokiol still enhanced the PGI<sub>2</sub> level in endothelial cells media in the presence of prostacyclin synthase inhibitor, but did not influence the effect of COX inhibitor, which suggested that the target of honokiol's effect was prostacyclin synthase not COX. The opinion that the target enzyme of honokiol's action was prostacyclin synthase seems disagreed with the idea that honokiol inhibited TXA<sub>2</sub> formation in platelets (Teng et al., 1988). One possible explanation is that honokiol inhibits TXA<sub>2</sub> synthesis in platelets, and augments PGI<sub>2</sub> generation in endothelial cells.

In addition, honokiol up-regulated prostacyclin synthase protein levels in endothelial cells in this study (Fig. 4). Honokiol




Fig. 6. Effect of honokiol on the caspase-3 activity. Rat aortic endothelial cells of passage 5 were seeded at  $2 \times 10^6$  cells/ml in 24-well plates and grown to confluence. After pre-incubated with vehicle or honokiol for 30 min, cultures were exposed to ox-LDL for 24 h to induce apoptosis, and in the other cultures, cells were only incubated with honokiol to evaluate honokiol per se effect on apoptosis of endothelial cells. The activity of caspase-3 was measured spectrophotometrically at 405 nm using a commercially available kit.  $OD_{405\text{treatment}}/OD_{405\text{control}}$  was used to express the increase in caspase-3 activity. Results were depicted as mean  $\pm$  S.D. \*\* $P < 0.01$ .

was certainly an enhancement factor for PGI<sub>2</sub> generation. Therefore the stimulation of honokiol on PGI<sub>2</sub> production is attributed to the increase of prostacyclin synthase protein expression.

Abnormalities of endothelial cell structure and function may contribute to the occurrence of diseases such as thrombosis and atherosclerosis. There are two mechanisms of cell death: apoptosis and necrosis. ox-LDL can induce apoptosis in endothelial cells (Lin et al., 2004; Zhang et al., 2003). Honokiol protected endothelial cells against ox-LDL injury (Hu et al., 2005) and showed significant anti-oxidation effect (Lo et al., 1994). To further reveal the mechanisms of its cell protection, the effect of honokiol on ox-LDL-induced apoptosis was explored in this study. The results of flow cytometry analysis and caspase-3 activity assay indicated the inhibitory effect of honokiol on the endothelial cell apoptosis was dependent on the concentration. 0.376–3.76 μM of honokiol significantly inhibited cell apoptosis, whereas a higher concentration (37.6 μM) failed to show this effect (Figs. 5 and 6). However, honokiol at all these concentrations (0.376–37.6 μM) did not change the activity of caspase-3 of normal endothelial cell (Fig. 6); suggesting that honokiol did not cause the apoptosis of normal cells at these concentrations.

It was reported that 40 μM of honokiol induced human squamous lung cancer CH27 cells apoptosis (Yang et al., 2002). This indicates that at high concentrations, honokiol may influence cells differently. In our previous study, it was found that the protective effect of honokiol at 37.6 μM on endothelial cells was not more potent than at 3.76 μM (Hu et al., 2005). These manifestations may be due to its maximum effect having been reached at 3.76 μM, indicating that at much higher concentrations (37.6 μM or higher) activity would not increase or decrease. There are various mechanisms of cell protection including inhibiting apoptosis, anti-oxidation and scavenging free radicals, etc. Based on the data of this study, the effect of honokiol was dependent on concentration, at lower concentration the inhibition of apoptosis was the central mechanism for protective effect of honokiol, while at higher concentrations the other factors were important (such as anti-oxidation). In accordance to the results it can be explained that 0.376–3.76 μM of honokiol significantly inhibited cell apoptosis, but higher concentrations (37.6 μM) failed. However, the detailed mechanism of honokiol on cell protection needs to be further explored.

In conclusion, honokiol augments PGI<sub>2</sub> generation in normal endothelial cells; its effect on PGI<sub>2</sub> generation attributes to up-regulation of prostacyclin synthase expression; its cell protection is correlated with its inhibition on apoptosis of endothelial cell. These findings have partly revealed the molecular mechanism of honokiol on inhibiting arterial thrombosis.

## Acknowledgments

Special thanks to Dr. Zhiyu Tang for the technical assistance. This study is supported by finance of National High-tech Research and Development Plan (863) No.2002AA2Z343C.

## References

- Centra, M., Ratych, R.E., Cao, G.L., Li, J., Williams, E., Taylor, R.M., Rosen, G.M., 1992. Culture of bovine pulmonary artery endothelial cells on Gelfoam blocks. *FASEB J.* 6, 3117–3121.
- Chaouchi, N., Wallon, C., Taib, J., Auffredou, M.T., Tertian, G., Lemoine, F. M., 1994. Interferon- $\alpha$ -mediated prevention of in vitro apoptosis of chronic lymphatic leukemia B cells: role of bel-2 and c-myc. *Clin. Immunol. Immunopathol.* 73, 197–204.
- Chen, S.Z., Wang, H., Wang, Y.Y., Hu, Y.H., 2003. Application of Honokiol, Magnolol and their preparation in cardiovascular and cerebral vascular diseases. Chinese patent: 200310121303.0.
- Chiu, J.H., Ho, C.T., Wei, Y.H., Lui, W.Y., Hong, C.Y., 1997. In vitro and in vivo protective effect of honokiol on rat liver from peroxidative injury. *Life Sci.* 61, 1961–1971.
- Cines, D.B., Pollak, E.S., Buck, C.A., Loscalzo, J., Zimmerman, G.A., McEver, R.P., 1998. Endothelial cells in physiology and in the pathophysiology of vascular disorders. *Blood* 91, 3527–3561.
- Debbie, B., Mary, E.M., 1995. Interleukin 1 induces prostacyclin dependent increases in cyclic AMP production and does not affect cyclic GMP production in human vascular smooth muscle cells. *Cytokine* 7, 417–426.
- Feng, L., Sun, W., Xia, Y., 1993. Cloning of two isoforms of rat cyclooxygenase: differential regulation of their expression. *Arch. Biochem. Biophys.* 307, 361–368.
- Golino, P., Ashton, J.H., Buja, L.M., Rosolowsky, M., Taylor, A.L., McNatt, J., Campbell, W.B., Willerson, J.T., 1989. Local platelet activation causes vasoconstriction of large epicardial canine coronary arteries in vivo. Thromboxane A<sub>2</sub> and serotonin are possible mediators. *Circulation* 79, 154–166.
- Hu, H., Zhang, X.X., Wang, Y.Y., Chen, S.Z., 2005. Honokiol inhibits arterial thrombosis through endothelial cell protection and stimulation of prostacyclin. *Acta Pharmacol. Sin.* 26, 1063–1068.
- Jaffe, E.A., Hoyer, L.W., Nachman, R.L., 1973. Synthesis of antihemophilic factor antigen by cultured human endothelial cells. *J. Clin. Invest.* 52, 2757–2764.
- Kishi, Y., Numano, F., 1989. In vitro study of vascular endothelial injury by activated platelets and its prevention. *Atherosclerosis* 76, 95–101.
- Lin, S.J., Shyue, S.K., Liu, P.L., Chen, Y.H., Ku, H.H., 2004. Adenovirus mediated overexpression of catalase attenuates ox-LDL induced apoptosis in human aortic endothelial cells via AP-1 and C-Jun N-terminal kinase/extracellular signal-regulated kinase pathway. *J. Mol. Cell. Cardiol.* 36, 129–139.
- Liou, K.T., Shen, Y.C., Chen, C.F., Tsao, C.M., Tsai, S.K., 2003. Honokiol protects rat brain from focal cerebral ischemia-reperfusion injury by inhibiting neutrophil infiltration and reactive oxygen species production. *Brain Res.* 992, 159–166.
- Lo, Y., Teng, C.M., Chen, C.F., Chen, C.C., Hong, C.Y., 1994. Magnolol and honokiol isolated from *Magnolia officinalis* protect rat heart mitochondria against lipid peroxidation. *Biochem. Pharmacol.* 47, 549–553.
- Ottino, Paulo, Duncan, John, R., 1997. Effect of vitamin E succinate on free radical formation, lipid peroxidation levels and cyclooxygenase activity in murine melanoma cells. *Nutr. Res.* 17, 661–676.
- Romualdo, J., Seguro, J., Babalola, O., Chieko, Y., Tadashi, T., 1997. Cyclic strain is a weak inducer of prostacyclin synthase expression in bovine endothelial cells. *J. Surg. Res.* 69, 135–318.
- Sibel, Ulker, Mehtap, G., 2001. Aprotinin impairs endothelium-dependent relaxation in rat aorta and inhibits nitric oxide release from rat coronary endothelial cells. *Cardiovasc. Res.* 50, 589–596.
- Spisni, E., Bartolini, G.H., Orlandi, M., 1995. Prostacyclin (PGI<sub>2</sub>) synthase is a constitutively expressed enzyme in human endothelial cells. *Exp. Cell Res.* 219, 507–513.
- Teng, C.M., Chen, C.C., Ko, F.N., Lee, L.G., Huang, T.F., Chen, Y.P., Hsu, H.Y., 1988. Two anti-platelet agents from *Magnolia officinalis*. *Thromb. Res.* 50, 757–765.
- Towbin, H., Staehlin, T., Gordon, J., 1979. Electrophoretic transfer of proteins from polyacrylamide gels to nitrocellulose sheets: procedure and some applications. *Proc. Natl. Acad. Sci. U. S. A.* 76, 4350.
- Watanabe, K., Watanabe, H., Goto, Y., Yamaguchi, M., Yamamoto, N., Hagino, K., 1983. Pharmacological properties of magnolol and honokiol extracted

- from *Magnolia officinalis*: central depressant effects. *Planta Med.* 49, 103–108.
- Watanabe, R., Kishi, Y., Sakita, S., Numano, F., 1997. Impaired NO release from bovine aortic endothelial cells exposed to activated platelets. *Atherosclerosis* 128, 19–26.
- Yang, S.E., Hsieh, M.T., Tsai, T.H., Hsu, S.L., 2002. Down-modulation of Bcl-XL, release of cytochrome *c* and sequential activation of caspases during honokiol-induced apoptosis in human squamous lung cancer CH27 cells. *Biochem. Pharmacol.* 63, 641–651.
- Zhang, W.G., Madhu, S., Li, D.Y., Jawahar, L., 2003. Metha. Identification of apoptosis-inducing factor in human coronary artery endothelial cells. *Biochem. Biophys. Res. Commun.* 301, 147–151.